EN

**DeVilbiss® SmartLink™ System  
DV5SS Software & DV5M Module**

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## TABLE OF CONTENTS

IMPORTANT SAFEGUARDS .....	3
INTRODUCTION	
Indications For Use .....	3
Compatibility .....	3
Overview .....	3
QUICK START GUIDE	
Software .....	4
Hardware .....	4
Data Transfer .....	4
SOFTWARE	
Installing & Uninstalling .....	5
Navigating .....	6
Software Tools .....	7
MANAGING THERAPY & DATA	
Managing Therapy Settings .....	9
Collecting Therapy Data .....	10
Storing Therapy Data .....	10
Transferring Data .....	10
Reporting Therapy Data .....	11
Generating Therapy Reports .....	12
MANAGING RECORDS	
Patient Records .....	13
Device Records .....	17
Study Records .....	20
TROUBLESHOOTING	
Software .....	23
Module .....	23
PARTS .....	24
SPECIFICATIONS .....	24
DECLARATION OF CONFORMITY .....	25
LICENSE .....	26
WARRANTY .....	26
APPENDIX	
Record Setup Screens .....	27
Therapy Setting Options and Defaults .....	29
Sample Reports .....	30
Patient Notes—System Notes (Auto-Entry) .....	34
SmartLink Messages (shown on CPAP LCD Display) .....	35

## IMPORTANT SAFEGUARDS

When using electrical products, basic safety precautions should always be followed. Important safety information in this manual is highlighted by the following terms:

- CAUTION** Information for preventing damage to the product.
- NOTE** Information to which you should pay special attention.

**PLEASE READ ALL INSTRUCTIONS BEFORE USING THIS DEVICE.**

**SAVE THESE INSTRUCTIONS.**

**CAUTION**—Do not use the DeVilbiss® SmartLink™ Desktop Software on computers with Operating Systems (OS) that have not been tested by DeVilbiss for compatibility. Only Microsoft® Windows® XP™ or Microsoft® Windows VISTA® operating systems are compatible. Using the software with an incompatible OS may cause the data to display incorrectly.

**CAUTION**—The circular data port connectors located on the back of the DeVilbiss SmartLink Module and the DeVilbiss IntelliPAP and SleepCube devices should only be used to attach the module to the CPAPs. Do not attempt to attach other devices to either the module or the CPAPs via these connectors as it may damage the module, the CPAP, or the device you are attempting to connect.

**CAUTION**—Use DeVilbiss SmartLink Data Cards with the SmartLink System to ensure correct data collection and display. Other Secure Digital (SD) cards may not be compatible.

## INTRODUCTION

### Indications For Use

The DeVilbiss SmartLink System can only be used in conjunction with the DV5I and DV54 Series CPAP Systems for follow up of obstructive sleep apnea patients weighing above 30 kg on nasal CPAP therapy.

### Compatibility

The SmartLink System is compatible with the following DeVilbiss CPAPs:

- IntelliPAP™ Standard (DV5I series) and IntelliPAP™ AutoAdjust® (DV54 series)
- SleepCube™ Standard (DV5I series) and SleepCube™ AutoAdjust® (DV54 series)

### Overview

The SmartLink System provides access to patient therapy information through the two-way transfer of data between patient devices and the SmartLink software including therapy efficacy and CPAP settings changes as requested by the attending physician.

The SmartLink System consists of the SmartLink Module, SmartLink Desktop Software, data card, and card reader.

- **SmartLink Desktop Software** is a computer program that stores and displays data such as patient, device, and study records. The software collects this data both from direct manual input and data transfers from the CPAP. The software also transfers therapy setting changes to the CPAP when needed.
- **SmartLink Module** is a communication device that attaches to the CPAP and collects therapy data. The module transfers information to and from the SmartLink software via a data card and card reader.
- **SmartLink Data Card** is a portable memory card that carries information from one data transfer device to another. In this system, the data card is a Secure Digital (SD) card that carries therapy information between the SmartLink module on the CPAP and the SmartLink software.
- **Card Reader** is a data transfer device connected to a computer that allows the transfer of data from a data card to the software and visa versa.
- **Oximetry** (optional) is a non-invasive measurement of the oxygen saturation levels of hemoglobin.

**CPAP TO Software** — Therapy data generated in the CPAP transfers to the SmartLink software.

CPAP w/SmartLink Module → Data Card → Card Reader → SmartLink Software

**Software TO CPAP** — Setting changes generated in the SmartLink software transfer to the CPAP.

SmartLink Software → Card Reader → Data Card → CPAP w/SmartLink Module

## QUICK START GUIDE

### Software

#### Installing software and drivers

1. Ensure you have Administration rights and privileges to the computer.
  2. Insert SmartLink CD in the computer's CD-ROM drive.
  3. Follow SmartLink Setup Wizard instructions.
  4. While logged in as Administrator, launch the application from the desktop icon.
  5. If prompted to allow the application to pass through the firewall, select ALLOW.
- All drivers are pre-installed in Microsoft® Windows XP® and Vista® operating systems.

#### Setting up study records

1. On the Patient table in the SmartLink software, click the 'New Patient' button.
2. Complete patient information and click SAVE.
3. When the Assign Device? prompt appears, select YES.
4. Select 'New Device' in the Select a Device menu.
5. Enter the serial number for the new device and click SAVE. A message will appear saying, Device created successfully.
6. In the Study Details screen, modify or accept the default settings in the Therapy Settings tab, the Device Settings tab and, for the AutoAdjust device only, the Respiratory Event Detection tab.
7. Click SAVE.
8. You will be prompted to transfer the settings to a data card, select 'OK' to proceed or 'Cancel' to verify that you will manually update device settings.

### Hardware

#### Installing the module

1. Align the module's primary data connector and its three locking tabs with the circular connector and the three locking recesses on the CPAP.
2. Press the module into the back of the CPAP until the module audibly clicks into a locked position.

#### Installing the data card in the module

1. Insert the data card into the slot on the top of the module until it is fully seated.

#### Installing the card reader and data card

1. Insert the USB connector on the card reader into the USB port on the computer.
2. Insert the card into the slot on the card reader.

### Data Transfer

#### Transferring from CPAP to software

1. Insert a data card into the slot on the top of the SmartLink™ Module until it is fully seated.
2. The CPAP will display one of the following messages while it attempts to transfer data:
  - Card Detected, Please Wait... indicates that the module is attempting to transfer data.
  - Card Detected, Ready for Use indicates that the card contains all current therapy data and is ready to receive new data as it becomes available.
  - Card Error, Contact Provider indicates that no data can be transferred.
3. After the transfer attempt, the CPAP will display one the following messages:
  - Card Transfer, Successful indicates that new data on the card transferred to the module and new data on the module transferred to the card. The data card can be removed or it can remain in the module to collect new therapy data as it becomes available from the CPAP. If you are removing the card from the module, press down on the data card to release it from the slot—the CPAP will display Card Removed, Send to Provider.
  - Card Error, Contact Provider indicates that the transfer of data to the card was not successful.
  - Setting Error, Contact Provider indicates that transfer of setting info was not successful. **NOTE—Provider intervention is necessary to confirm prescription settings are correct.**
4. Attach a card reader to the computer running SmartLink software.
5. Insert the data card holding the new therapy data into the card reader.
6. One of two prompts will appear in the software: Create a new study...? or Upload new data? Follow software instructions for your chosen prompt. Data will download automatically.

## SOFTWARE

**IMPORTANT: ONLY a recognized Power User or Administrator in Microsoft® Windows XP® or Vista® operating systems can run and use the DeVilbiss SmartLink Desktop Software.**

### Installing & Uninstalling

#### System Requirements

Operating System	Microsoft® Windows® XP™ service pack 2 or higher OR Microsoft® Windows® Vista™ service pack 1 or higher
Adobe® Reader®	6.0 or greater
RAM	512 MB or greater
Hard Drive Free Space	150 MB
Hardware:	SD Memory card reader supporting FAT16 file format OR USB port for auxiliary card reader

#### SmartLink Desktop Software

This software is not intended for use on a network. Only a single instance of the application can access the database at one time.

Software files are shipped in a compressed format and will install on your hard drive in directory C:\Program Files\SmartLink Desktop unless you specify otherwise when prompted by the setup program. You must install the software, not merely copy it to the hard drive.

#### Installing SmartLink software

1. Make sure you have Administrator privileges and write permission in the following locations:

- Folder where you are installing SmartLink Software
- Microsoft® Windows system folder
- System registry

You must have permission to write to the necessary system registry sections. Refer to Using Local Users and Groups in your system's online Help, which you can access by clicking Start on the taskbar.

2. Ensure that you meet the minimum system requirements listed above.
3. Disable all anti-virus programs operating on your system and close any open applications.
4. Insert the installation disc into your CD-ROM drive.
5. SmartLink™ Setup Wizard will start automatically. Follow on-screen instructions.

**NOTE**—If the wizard does not start automatically: Open the file: d:\setup.msi (where 'd' is the CD-ROM drive letter) to force the wizard to start.

6. While logged in as Administrator, launch the application from the desktop icon.
7. If prompted to allow the application to pass through the firewall, select ALLOW.

#### Uninstalling SmartLink software

Uninstalling removes all files that were installed using the SmartLink CD. Files created after SmartLink was installed, such as database and exported data, are not removed. You can manually delete these files after uninstalling the software.

1. From the Windows Start menu, select 'All Programs'.
2. Select 'SmartLink Desktop'.
3. Select 'Uninstall'.
4. When the uninstall is complete, restart your computer to ensure that any files that may have been in use by your operating system are properly uninstalled. Select 'Cancel' if you do not want to restart your computer immediately.

#### Card Reader Software and Drivers


The card reader requires no software. It is a plug-and-play accessory. All drivers are pre-loaded in the Microsoft® Windows® XP™ and Vista™ operating systems.

- If more than one SmartLink data card is installed in the computer, a warning appears: Multiple removable media devices are attached. Disconnect or eject all of the following drives to continue: x:\ y:\ (where x:\ and y:\ are the extra data cards).


## Navigating


Five menus are available on the main SmartLink software page: Languages, Patients, Devices, Administration, and My Profile.


**Languages menu** contains five available options: 'English', 'French', 'German', 'Spanish', and 'Italian'.


**Patients menu** contains five tables holding patient information: 'All Patients', 'New Studies', 'Studies in Process', 'Studies Complete', 'Archived Patients'. More information on these tables is available in [Software / Software Tools / Working with Tables](#).

**Devices menu** contains three tables holding device information: 'All Devices', 'Available Devices', 'Assigned Devices'. More information on these tables is available in [Software / Software Tools / Working with Tables](#).

**Administration menu** contains a link to the Patient Access Log. More information on this log is available in [Software / Software Tools / Patient Access Log](#).

**My Profile** contains a link to the User Information form. More information on this form is available in [Software / Software Tools / My Profile](#).

## Software Tools

### Working with Tables

The tables found in the Patients and Devices menus contain tools to help you navigate through large quantities of records.

SmartLink Desktop

Languages

SmartLink

DeVilbiss Healthcare

Patients Devices Administration System Administrator [My Profile]

All Patients

1 New Patient 2 10 items per page

4 19 items found, displaying 1 to 10. [First/Prev] 1, 2 [Next/Last] 5

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	« Select an option: »
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	« Select an option: »
Rlnnm	Rnhert	<Unassigned>					« Select an option: »
Casey	Kathy	IntelliPAP AutoAdjust	HD003488			New Study	« Select an option: »
Collins	Mike	<Unassigned>					« Select an option: »
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	« Select an option: »
Dell	Jamie	<Unassigned>					« Select an option: »
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	« Select an option: »
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	« Select an option: »
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	« Select an option: »

Version 0.1.40 © 2008 DeVilbiss Healthcare

1. **[New Patient] and [New Device] buttons** allow you to add items to the table.
2. **Page display option** allows you to select how many items to show per page. Possible values are 10, 25, 50, All.
3. **Sort selector** allows you to sort columns by clicking in the column header. Reclicking a column will reverse the order.
4. **Inventory display** shows the total number of records available in this table and the number of records displayed on this page of the table.
5. **Page selector** allows you to move to the first and last pages, to previous and next pages, or to a specific page number in multi-paged table.
6. **Options menu** takes you directly to specific record views or action windows.

## Patient Access Log

SmartLink software tracks access to patient records in the Patient Access Log. Every day that study or patient records are opened, an entry is added to the access log. Only one entry per patient is added on any one day, even if a user opens the same patient's record multiple times during that day.

1. Select 'Patient Access Log' under the Administration menu.

Patients	Devices	Administration	System Administrator [My Profile]
All Patients		Patient Access Log	

2. A table of log entries appears.

**NOTE:** Click on column heading to sort the entries.

Patients	Devices	Administration	System Administrator [My Profile]
Patient Access Log			
<div>10 items per page</div> <div>3 items found, displaying all items.</div> <div>1</div>			
Date	Patient	User Name	User Email
7/10/2008	R. U. Sleepy	System Administrator	
7/10/2008	Rip Van Winkle	System Administrator	
7/10/2008	Little Winken	System Administrator	

## My Profile

SmartLink software allows you to enter information unique to your company in a User Profile form.

1. Enter or modify company information as you want it to appear on all printed and exported reports.
2. Enter 'Low Use Threshold' hours—the minimum amount of time the patient must breathe on the CPAP to be counted as a full night. This entry will determine how usage is computed.

Required Fields are marked

User Information	
Name	System Administrator
Address	
City	
State/Province	
Postal Code	
Country	
Phone Number	
Email Address	
Contact	
Low Use Threshold	4 Hours
Percentile Pressure	90 Percentile

Save Cancel

3. Enter 90th or 95th Percentile Pressure – The selected pressure will be displayed on the appropriate reports.
4. Click SAVE.

## MANAGING THERAPY AND DATA

### Managing Therapy Settings

DeVilbiss® DV5X series CPAPs and SmartLink™ Therapy Management provide you with a wide variety of therapy setting options. Full details on the options and their default settings are listed in the [Appendix – Therapy Setting Options and Defaults](#). You can set therapy parameters manually on the CPAP itself or set them in the SmartLink software and then transfer them to the CPAP using a data card.

To manually set therapy parameters on the CPAP, follow the instructions in the DV51 or DV54 Set-Up Guides. Verify that the settings match those in the SmartLink software study record.

**To set therapy parameters in SmartLink software and then transfer the new settings to the CPAP using a data card:**

1. Insert the card reader into a USB port on the computer running SmartLink software and insert a data card in the reader.
2. Select the appropriate table in the Patients menu.
3. Select 'View/Edit Study Details' from the Options menu for that patient.

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: >
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	View Reports Export Oximetry Data Update Patient Notes View/Edit Study Details View/Edit Patient Information Un-assign Device Archive Patient
Dell	Jamie	<Unassigned>					< Select an option: >
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: >
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: >
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: >

4. Enter the desired changes in the Therapy Settings tab, the Device Settings tab or, for the AutoAdjust device only, the Respiratory Event Detection tab. Views of these study tabs and their default settings are available in the [Appendix](#).

5. Click the SAVE button to transfer the changes to the data card and to save them in the software.
  - If the transfer is successful, the software will enter a System Note in the patient record: SmartLink data card has been programmed with updated settings. Verification of settings is pending.
6. Insert the data card into the SmartLink Module attached to the patient's CPAP. The CPAP will display one of the following messages:
  - Card Detected, Please Wait... indicates that the module is attempting to transfer data. In addition to the setting changes sent from the software, the module itself may contain new therapy management data that it needs to transfer to the data card.
  - Card Transfer, Successful indicates that new data on the card transferred to the module and new data on the module transferred to the card. The data card can be removed or it can remain in the module to collect new therapy data as it becomes available from the CPAP.

**NOTE**—If the module contained new therapy study information, that data was transferred to the card after the new settings were transferred to the module and CPAP.

- Setting Error, Contact Provider indicates that transfer of setting info was not successful. **NOTE**—*Provider intervention is necessary to confirm prescription settings are correct.*
- Card Detected, Ready for Use indicates that the card contains all current therapy data and is ready to receive new data as it becomes available.
- Card Error, Contact Provider indicates that the transfer of data to the card was not successful.

### To set therapy parameters in SmartLink software and then setup the CPAP manually:

**NOTE**—*This action requires you to manually change the settings on CPAP.*

1. Follow steps 1 through 4 above.
2. Click SAVE
3. Click 'Cancel', when prompted to insert data card.
4. Click YES, when prompted to manually update settings on the CPAP.
5. Manually set therapy parameters on the CPAP.

## Collecting Therapy Data

When a SmartLink module is connected to an IntelliPAP or SleepCube, therapy data is automatically collected and stored in the module. This data includes information such as CPAP usage and oximetry if available. If the module is connected to an AutoAdjust model, it also collects information on CPAP pressures, leak, AHI, and respiratory events. Data collected in the module is transferred to the SmartLink software via a data card.

When you create a Patient Record, Device Record, or Study Record in the SmartLink software, the data is stored in the database. The data includes information that identifies the patient, device, and particular study parameters. Data transferred from the SmartLink module via a data card is added to this database.

## Storing Therapy Data

The SmartLink module has a rolling buffer capable of storing three years of usage hours, 120 nights of detailed performance, and 72 hours of oximetry data. When a data card is installed in the module, the module copies its information to the data card. The module does not erase data after a transfer to the data card; instead it stores all information collected from the CPAP until its buffer is full, then it begins to drop its oldest data to make space for new data.

The software stores information in a database. Information that is or has ever been associated with a study record cannot be deleted from the database. Patient and device records that were never assigned to study records can be deleted from the software database.

## Transferring Data

### Data Transfer from the CPAP to a Data Card

When the SmartLink module detects a data card, the module completes the following actions:

- The module scans the card for any changes to the therapy settings that might have been sent from the software. If changes are present, the module transfers new settings to the CPAP.
- The module checks its stored data against any data that may be on the card from previous transfers. If the module has new data, it is copied to the data card and the CPAP notifies the user that the transfer was successful.
  - The data card can be removed from the module after this transfer and sent to the clinician/provider for transfer to the patient's study record OR
  - The data card can be kept in the module and removed at a later date.

### Data Transfer from the Data Card to the Software

When the software detects a data card in an attached card reader, the software prompts you to complete the following actions:

- If the data card contains data from a CPAP that is not recognized as part of a study record, the Create a new study? prompt appears. Select YES to complete the new study setup and transfer the data from the card to the study record.
- If the data card contains data from a CPAP that is recognized as part of a study record, the Upload new data? prompt appears. Select YES to transfer the data from the card to the record.

## Reporting Therapy Data

**NOTE**—Oximetry reports are available separately through Nonin® nVISION® software. Use the instructions available in [Managing Therapy Data / Generating Therapy Reports](#).

There are two report styles available in SmartLink software for each patient's study record: Summary and One-Night.

Exported reports begin with a cover page that includes detailed patient information and a usage summary of the entire date range. The usage summary includes:

- Possible usage days in range
- Days with any usage
- Days used at least x hours (where x is equal to the low use threshold entered in My Profile)
- Usage index (days used at the low use threshold or higher/possible usage days)
- Average usage duration
- Standard deviation

## One-Night Reports


The One-Night report shows detailed therapy data for each night of recorded use. A night is defined as at least 30 minutes of continuous use. Use the Sessions menu to view the report for the chosen night.

**NOTE**—If an oximeter was used while sleep therapy was delivered by a DV54 AutoAdjust CPAP, the One-Night report will show that activity. To create a full report in Nonin® nVISION® software, use the instructions available in [Managing Therapy Data / Generating Therapy Reports](#).

## Summary Reports

The Summary report shows averaged therapy data for a specific date range. The Summary report is available in two views: Summary and Detail.

In Summary view, the usage is shown as the total number of hours the patient used the CPAP during a 24 hour period (noon to noon) in the selected range of usage.

In the Detail view, the usage is shown as the ON/OFF times for each 24 hour period (noon to noon) in the selected range of usage indicating when during each day and for how long the patient used the CPAP.

In both views of the Summary reports, a specific range of nights can be requested as long as it is within the date range for the entire body of study data collected.

## Generating Therapy Reports

**NOTE**—Reports are only available after therapy data has been transferred to the software. Use the instructions available in [Managing Therapy Data / Transferring Data](#) to transfer therapy data to the software.

**To generate a SmartLink™ Therapy Management report, follow these steps:**

1. Click the Patients menu and select 'Studies in Process' or 'Studies Complete' to locate the chosen patient for the report.
2. Select 'View Reports' in the Options menu for the chosen patient.
3. Click either the 'Summary' or 'I-Night' tab, to select the report format of your choice. The report will refresh after this change.

**NOTE**—A Summary report has two additional options: the Detail view contains a graph showing ON/OFF usage for each night and the Summary view contains a graph showing the total hours of usage for one or more nights.

4. Select the date or date range for the report.

**NOTE**—The dates and date ranges must be within the dates of the transferred therapy data.

- For Summary reports, enter the report date range in the 'From' and 'To' boxes. The report will refresh after this change.
  - For I-Night reports, open the Sessions menu and select the night for the report. The report will refresh after this change.
5. Optional: If printing or sending the report to another person or file, click the [Export to PDF] button; then print the report from the PDF screen or save it to a file for electronic transfer.

**To generate an Oximetry report, follow these steps:**

1. Ensure that Nonin® nVISION® software is installed on your system.
2. Open SmartLink™ Desktop Software.
3. Click the Patients menu and select 'Studies in Process' or 'Studies Complete' to locate the chosen patient for the report.
4. Select 'Export Oximetry data' in the Options menu for the chosen patient. A message will appear saying Loading into sessions.
5. When the Software Installation Path prompt appears, select OK. The default storage location is C:\Program Files\SmartLink Desktop\Nonin.
6. Open nVISION® software and follow the software guide to open the files and generate reports.

## Sample Reports

Views of sample reports for the DV54 AutoAdjust CPAP and the DV5I Standard CPAP are available in the [Appendix](#).

## MANAGING RECORDS

### Patient Records

#### Creating a Patient Record

Patient records hold a variety of useful information about the patient such as contact, payor, and physician information. You can create a patient record as part of a new study or as a stand alone record that can be used in a study later.

1. Select 'All Patients' in the Patients menu; then click the [New Patient] button.

The screenshot shows a software interface with a top navigation bar containing 'Patients', 'Devices', 'Administration', and 'System Administrator [My Profile]'. A dropdown menu is open under 'Patients', showing options: 'All Patients', 'New Studies', 'Studies in Process', 'Studies Complete', and 'Archived Patients'. The 'All Patients' option is highlighted. In the main content area, there is a 'New Patient' button and the text 'Nothing found to display.'

2. Complete patient information and click SAVE. A view of the patient information form is available in the [Appendix](#).

**NOTE:** You can create a record using only the patient's first and last name. Additional information can be added to the record at any time.

3. When the Assign Device? prompt appears,
  - If you are creating a stand-alone record, click NO to finish.
  - If you are setting up a new study, select YES to proceed; continue with Step 4.
4. In the first Study Details window, click the 'Select a Device' menu and continue with one of the following options:
  - Select [New Device] and enter the serial number for the new device OR
  - Select the serial number of the device to be used in the study

The screenshot shows the 'Study Details' window. It has a title bar 'Study Details'. Inside, there are four sections: 'Select a Patient' with a dropdown menu showing 'Snore NoMore' and a note '\* Archived patients will be restored upon save.'; 'Select a Device' with a dropdown menu showing 'Select a Device'; 'Study Length' with a dropdown menu showing 'Select a Device' and a text input field; and 'Comments' with a text area. At the bottom, there are 'Save' and 'Cancel' buttons.

5. Click SAVE. A message will appear saying, Device created successfully.
6. In the next Study Details screen, modify or accept the default settings in the Therapy Settings tab, the Device Settings tab and, for the AutoAdjust device only, the Respiratory Event Detection tab. Views of these study tabs and their default settings are available in the [Appendix](#).
7. Click SAVE.
8. When prompted to insert a data card to transfer the study settings to the CPAP:
  - Select Cancel to configure the CPAP settings manually.
  - Select OK to copy the study settings to the data card for transfer to the CPAP through the SmartLink Module.

**NOTE**—A study's status is listed as New Study until the first therapy data transfers from the CPAP to the software; then the status changes to Study in Process. More information on study status is available in [Managing Records / Study Records](#).

## Patient Notes

Patient Notes hold valuable information about a patient's therapy including data transfers, events, and observations. All notes appear on therapy reports and contain the date and source of entry. Patient Notes are added automatically by the management system (System Notes) and can be added manually by the clinician/provider (User Notes). System Notes cannot be edited or deleted.

### System Notes are entered when one of the following actions occurs:

- Successful data transfers are noted with one of the following:
  - First data received for study.
  - SmartLink data card has been programmed with updated settings. Verification of settings is pending.
  - Settings update occurred on <date of transfer to CPAP> and validated on <download date>.
  - New data has been uploaded.
- Errors in data transfer are noted with one of the following:
  - Download contains one or more faulty records.
  - Errors in data import. Some data may have been lost.
- Device assignments and un-assignments are noted with:
  - Device <serial number> was assigned to patient.
  - Device <serial number> was un-assigned from patient.

### User Notes are entered for information gathered during home visits or office calls. To create, edit or delete User Notes:

1. Select the appropriate table in the 'Patients' menu.
2. Find the chosen patient in the table and select 'Update Patient Notes' in the Options menu.

All Patients

New Patient

10 items per page

19 items found, displaying 1 to 10.  
[First/Prev] 1, 2 [Next/Last]

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: >
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: >
Dell	Jamie	<Unassigned>					< Select an option: >
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: >
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: >
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: >

3. Add, edit or delete notes:

- To ADD a note: click the [Add Patient Note] button, enter the note in the text box and click SAVE. This note will appear in the table and in patient reports.

Update Patient Notes - Kathy Casey

Add Patient Note

10 items per page

8 items found, displaying all items.

1

Date	Added By	Note	Options
8/20/08	System	First data received for study	N/A
8/20/08	System	Device HD003601 was assigned to patient	N/A

- To EDIT a note: select 'Edit' in the Options menu for the chosen note, edit the note in the text box and click SAVE.

Add Patient Note50 items per page

Enter Note Text Below:

Date7/11/08

Patient requested cables to connect PAP to battery for camping.

SaveCancel

- To DELETE a note: select 'Delete' in the Options menu for the chosen note and then click YES when prompted Are you sure you want to delete this patient note?

Viewing/Editing a Patient Record

- Select the appropriate table in the Patients menu.
- Find the chosen patient in the table and select 'View/Edit Patient Information' in the Options menu.

All Patients

New Patient10 items per page

19 items found, displaying 1 to 10.  
[First/Prev] 1, 2[Next/Last]

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: > View Reports Export Oximetry Data Update Patient Notes View/Edit Study Details View/Edit Patient Information Un-assign Device Archive Patient
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: >
Dell	Jamie	<Unassigned>					< Select an option: >
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: >
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: >
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: >

- If you are editing or adding information, click SAVE when you have finished changing the record. A view of the patient information form is available in the [Appendix](#).

Archiving a Patient Record

You may archive patient records that no longer require active therapy management. Archiving removes the record from the normal workflow. The record remains in the database, but is only visible through the 'Archived Patients' table in the Patients menu.

- Select 'All Patients' in the Patients menu.

- Find the chosen patient in the table and select 'Archive Patient' in the Options menu.

**All Patients**

[New Patient](#) 10 items per page

19 items found, displaying 1 to 10.  
[First/Prev] 1, 2[Next/Last]

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: >
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: >
Dell	Jamie	<Unassigned>					< Select an option: >
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: >
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: >
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: >

- When the Archive Patient? prompt appears, select YES to move the record to the 'Archived Patients' table. The record appears in the new table and no longer appears in any other table in the Patients menu.

### Restoring an Archived Patient Record

- Select 'Archived Patients' table in the Patients menu.

Patients	Devices	Administration	System Administrator [My Profile]
<a href="#">All Patients</a> <a href="#">New Studies</a> <a href="#">Studies in Process</a> <a href="#">Studies Complete</a> <a href="#">Archived Patients</a>		<a href="#">New Patient</a> Nothing found to display.	

- Find the chosen patient in the table and click [<<Restore Patient>>].

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Van Winkle	Rip	IntelliPAP AutoAdjust	HD001000			New Study	<< Restore Patient >>
Head	Sleepy	IntelliPAP AutoAdjust	HD001030			New Study	<< Restore Patient >>
Land	LaLa	IntelliPAP AutoAdjust	HD001070			New Study	<< Restore Patient >>

- The patient record moves from the 'Archived Patients' table and reappears in the appropriate tables.

Device Records


Creating a Device Record

Device records hold serial numbers and model types for CPAPs. You can create a device record as part of a new study or as a stand alone record that can be assigned to a study later.

To create a device record as part of a new study, use the instructions available in [Managing Records / Creating a Study Record](#).

To create a stand-alone device record:

- 1. In any table under the Devices menu, click the [New Device] button.

- 2. Enter the serial number and click SAVE. You can add comments to this record at any time.

A screenshot of the 'Device Details' form. It has a title bar 'Device Details' and a subtitle 'Required Fields are marked'. There are two main input areas: 'Serial Number' with a text box and a red asterisk, and 'Comments' with a larger text area. At the bottom are 'Save' and 'Cancel' buttons.

Editing a Device Record

Serial numbers in the device records cannot be edited; however, Comments can be added, deleted and modified.

- 1. Select the appropriate table under the Devices menu and find the chosen device.
- 2. Click [**<<Edit>>**] OR select 'Edit' in the Options menu.

A screenshot of the 'All Devices' table. The table has columns: Type, Name, Serial Number, Status, and Options. It lists 10 devices. The first device is 'DV510 IntelliPAP Standard' with serial number 'ID002500' and status 'Assigned to Jack Rhoades'. The last device is 'DV540 IntelliPAP AutoAdjust' with serial number 'HD00SP33' and status 'Available (Never Assigned)'. There is a 'New Device' button and a '10 items found, displaying all items.' message above the table. The 'Options' column has a dropdown menu open for the last device, showing 'Edit' and 'Delete' options.

- 3. Add, modify and delete text in the Comments box and click SAVE.

A screenshot of the 'Device Details' form. The 'Serial Number' field is filled with 'HD001000'. The 'Comments' field is empty. At the bottom are 'Save' and 'Cancel' buttons.

## Assigning a Device to a Patient

You can assign a device to a patient while you are creating a patient record or a device can be assigned to a previously created patient record. Assigning a device to a patient creates a new study.

To assign a device to a new patient, use the instructions available in [Managing Records / Creating a Patient Record](#) or [Creating a Study Record](#).

### To assign a device to an existing patient:

1. Select 'All Patients' in the Patients menu.
2. Find the chosen patient in the table and select 'Assign Device' in the Options menu.

All Patients

New Patient

10 items per page ▾

19 items found, displaying 1 to 10.  
[First/Prev] 1, 2[Next/Last]

Last Name ▾	First Name ▾	Device ▾	Serial Number	Start Date ▾	End Date ▾	Status ▾	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: > ▾
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: > ▾
Bloom	Robert	<Unassigned>					< Select an option: > ▾
Casey	Kathy	IntelliPAP AutoAdjust	HD003488			New Study	< Select an option: > ▾
Collins	Mike	<Unassigned>					< Select an option: > ▾
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: > ▾
Dell	Jamie	<Unassigned>					< Select an option: > ▾
Farmer	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: > ▾ Update Patient Notes View/Edit Patient Information Assign Device Archive Patient
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: > ▾
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: > ▾

3. Click the 'Select a Device' menu and select the serial number of the chosen device OR select [New Device].

Study Details	
Select a Patient	Little Blinken * Archived patients will be restored upon save.
Select a Device	Select a Device
Study Length	<div> <div>Select a Device</div> <div>[New Device]</div> <div>HD001001 (IntelliPAP AutoAdjust)</div> <div>HD001020 (IntelliPAP AutoAdjust)</div> </div>
Comments	
<div>Save</div> <div>Cancel</div>	

- a. If you selected a serial number, go to the next step.
  - b. If you selected [New Device], use the instructions available in [Managing Records / Creating a Study Record](#).
4. When the New Study window appears, modify or accept the default settings in the Therapy Settings tab, Device Settings tab and, for the AutoAdjust device only, the Respiratory Event Detection tab. Views of these study tabs and their default settings are available in the [Appendix](#).
  5. Select the number of days in the Study Length drop-down menu in the Therapy Settings tab.
  6. Enter pertinent comments now or add comments to this study later.
  7. Click SAVE

Un-Assigning a Device Record

To change a device during a study, the current device record must be un-assigned and a new device record assigned to the patient which will start a new study.

To un-assign a device record, follow these steps:

- 1. Select the appropriate table in the Patients menu.
- 2. Find the chosen patient in the table and select 'Un-Assign Device' in the Options menu.

All Patients

New Patient10 items per page19 items found, displaying 1 to 10.[First/Prev] 1, 2[Next/Last]

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: > View Reports Export Oximetry Data Update Patient Notes View/Edit Study Details View/Edit Patient Information Un-assign Device Archive Patient
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: >
Dell	Jamie	<Unassigned>					< Select an option: >
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: >
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: >
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: >

- 3. When the Un-Assign Device? prompt appears, click YES.

Deleting a Device Record

Devices that have never been assigned to a patient can be deleted from the software database. Devices that have been associated with study records, whether or not they collected and transferred therapy data, cannot be deleted from the database.

To replace a device assigned to an active study, use the instructions available in Managing Records / Un-assigning a Device Record and Assigning a Device to a Patient Record.

To delete a device record from the software, follow these steps:

- 1. Select 'Available Devices' in the Device menu.
- NOTE**–Only devices in the 'Available (Never Assigned)' status can be deleted.
- 2. Find the chosen device in the list and select 'Delete' from the Options menu.

Type	Name	Serial Number	Status	Options
DV54D	IntelliPAP AutoAdjust	HD001001	Available	<< Edit >>
DV54D	IntelliPAP AutoAdjust	HD001020	Available (Never Assigned)	< Select an option: > < Select an option: > Edit Delete

- 3. Select YES, when the Delete this device? prompt appears. A message will appear saying Device Deleted Successfully.

## Study Records

Study records hold patient records, device records, and pertinent study information such as therapy parameters and results. Therapy parameters can be added or modified in the study record at any time. A study begins by assigning a device to a patient.

Study records are assigned a status as follows:

- **New Study** status is assigned to all studies until they receive their first successful transfer of therapy data from the CPAP.
- **Study in Process** status is assigned to studies after the first successful data transfer. The study will remain in this status until the software receives a data transfer after the 'Length of Study' expired.
- **Study Complete** status is assigned to studies when therapy data is received with dates indicating that the information was collected after the 'Length of Study' had expired.

**NOTE**—The study record will continue to accept transfers of therapy data while in Study Complete status.

### Creating a Study Record Using EXISTING Patients and Devices

1. Select 'All Patients' in the Patients menu and find the chosen patient in the table.

**NOTE**—You can click the column headings to sort the table. This may allow you to find specific patients more easily.

2. Select 'Assign Device' in the Options menu.

All Patients

New Patient

10 items per page

19 items found, displaying 1 to 10.

First/Prev

1, 2Next/Last

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003488			New Study	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: >
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: >
Dell	Jamie	<Unassigned>					< Select an option: >
Farner	Glen	IntelliPAP AutoAdjust	HD001499			New Study	< Select an option: > Update Patient Notes View/Edit Patient Information Assign Device Archive Patient
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	< Select an option: >
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	< Select an option: >

3. In the first Study Details window, click the 'Select a Device' drop-down menu and either:
  - a. Select [New Device] and enter the serial number for the new device OR
  - b. Select the serial number of an existing device

Study Details	
Select a Patient	Jamie Dell * Archived patients will be restored upon save.
Select a Device	Select a Device
Study Length	[New Device] HD00SP33 (IntelliPAP AutoAdjust)
Comments	
<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

4. In the next Study Details screen, modify or accept the default settings in the Therapy Settings tab, the Device Settings tab and, for the AutoAdjust device only, the Respiratory Event Detection tab. Views of these study tabs and their default settings are available in the [Appendix](#).
5. Click SAVE.
6. When prompted to insert a data card to transfer settings to the CPAP:
  - a. Select OK to copy the settings to the data card for transfer to the CPAP through the SmartLink Module. – OR –
  - b. Select 'Cancel' and configure the CPAP settings manually.

## Creating a Study Record Using NEW Patients and Devices

1. Click the [New Patient] button OR select 'All Patients' in the Patients menu and then click the [New Patient] button.
  2. Enter the patient information and click SAVE. A view of the patient information form is available in the [Appendix](#).
- NOTE**—You can create a record using only the patient's first and last name. Additional information can be added to the record at any time.
3. When the Assign Device? prompt appears, select YES.
  4. In the Study Details window, select [New Device] in the 'Select a Device' menu and enter the device serial number in the next window.

5. Click SAVE. A message will appear saying Device Created Successfully.
6. In the next Study Details screen, modify or accept the default settings in the Therapy Settings tab, the Device Settings tab and, for the AutoAdjust device only, the Respiratory Event Detection tab. Views of these study tabs and their default settings are available in the [Appendix](#).

7. Click SAVE.
8. When prompted to insert a data card to transfer settings to the CPAP:
  - a. Select OK to copy the settings to the data card for transfer to the CPAP through the SmartLink Module. – OR –
  - b. Select 'Cancel' and configure the CPAP settings manually.

## Editing a Study Record

To change therapy settings in a study record and transfer the changes to the CPAP, use the instructions available in [Managing Therapy Data / Managing Therapy Settings](#).

To edit or add comments in a device record or to delete a device record, use the instructions available in [Managing Records / Editing a Device Record](#) or [Deleting a Device Record](#). You cannot make changes to device information through the study record.

### To edit patient or therapy information:

1. Select the appropriate table in the 'Patients' menu.
2. Click the Options drop-down menu for the chosen patient:
  - a. Select 'View/Edit Study Details' to edit therapy information.
  - b. Select 'View/Edit Patient Information' to edit patient details.

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	< Select an option: >
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	< Select an option: >
Bloom	Robert	<Unassigned>					< Select an option: >
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	< Select an option: >
Collins	Mike	<Unassigned>					< Select an option: >
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	< Select an option: >
Dell	Jamie	<Unassigned>					< Select an option: >

**NOTE**—If you selected 'View/Edit Study Details', use the tabs under [Device Settings] button to locate the information you want to edit in the first Study Details screen. Views of these tabs and their default settings are available in the [Appendix](#).

- In the next Study Details screen, edit the settings in the Therapy Settings tab, the Device Settings tab and, for the AutoAdjust device only, the Respiratory Event Detection tab. Views of these study tabs and their default settings are available in the [Appendix](#).

- Click SAVE.
- When prompted, insert a data card to transfer settings to the CPAP:
  - Select OK to copy the settings to the data card for transfer to the CPAP through the SmartLink Module. – OR –
  - Select 'Cancel' and configure the CPAP settings manually.

### Working with Study Records

- Select the appropriate table in the 'Patients' menu.
- Click the Options menu for the chosen patient.

**NOTE**—Only those options which are appropriate for the current study status are available.

All Patients

New Patient

10 items per page

19 items found, displaying 1 to 10.  
[First/Prev] 1, 2 [Next/Last]

Last Name	First Name	Device	Serial Number	Start Date	End Date	Status	Options
Atkins	Judy	IntelliPAP AutoAdjust	HD003226			New Study	« Select an option: »
Blair	Shawn	IntelliPAP Standard	ID001444			New Study	« Select an option: »
Bloom	Robert	<Unassigned>					« Select an option: »
Casey	Kathy	IntelliPAP AutoAdjust	HD003601	7/25/08	9/23/08	Study in Process	« Select an option: »
Collins	Mike	<Unassigned>					« Select an option: » View Reports Export Oximetry Data Update Patient Notes View/Edit Study Details View/Edit Patient Information Un-assign Device Archive Patient
Conner	David	IntelliPAP AutoAdjust	HD005423			New Study	« Select an option: »
Dell	Jamie	<Unassigned>					« Select an option: »
Farrer	Glen	IntelliPAP AutoAdjust	HD001499			New Study	« Select an option: »
Fuller	Mary	IntelliPAP Standard	ID001554			New Study	« Select an option: »
Gibbs	Donald	IntelliPAP Standard	ID003598			New Study	« Select an option: »

- View Reports** – Links directly to views of graphs, charts and other data displays derived from a specific patient's therapy data. These views can be saved and printed as study reports in PDF format. The View Reports option is only available after therapy data is downloaded.
- Export Oximetry Data** – Initiates copying Oximetry data into session files that will default to storage in C:\Program Files\SmartLink Desktop\Nonin. This option is only available if oximetry data was transferred to the SmartLink software for this patient's study record. You must use Nonin® nVISION® software to view the data.
- Update Patient Notes** – Links directly to the list of System and User Notes for the chosen patient. This option allows you to add, delete or edit user notes. You cannot edit or delete system notes from patient records.
- View/Edit Study Details** – Links directly to the study record and allows you to change therapy settings. This option is only available after creating a study.
- View/Edit Patient Information** – Links directly to the patient record and allows you to add or change information concerning the patient.
- Assign Device** – Allows you to select a current device record or create a new device record that will be automatically assigned to the patient. When you assign a device to a patient, you create a study record.
- Un-Assign Device** – Links directly to a prompt asking if you want to un-assign the device. Selecting Yes removes the assigned device from the study record, terminates the study. If a new device is assigned, a new study begins.
- Archive Patient** – Links directly to a prompt asking if you want to archive this patient. Selecting Yes removes the patient record from all tables in the Patient menu except 'Archived Patients'. Archived patients can be restored. Patients that do not require further data management should be archived.

## TROUBLESHOOTING

The DeVilbiss SmartLink Module has no serviceable parts.

## Software

Problem	Probable Cause	Remedy
Software Installation problems	Computer does not meet minimum requirements.	Make sure computer meets minimum system requirements of memory, hard disk space, and operating system listed in Software Installing and Uninstalling section.
Software does not start	Windows Vista firewall blocking SmartLink software from running.	After installing on Windows Vista, SmartLink needs to be opened once by the Administrator to allow it to run.
	Computer does not meet minimum requirements.	Make sure computer meets minimum system requirements of memory, hard disk space, and operating system listed in Software Installing and Uninstalling section.
Software locks up	Not enough system resources to run SmartLink software.	Close other Windows programs to free system resources.
	Data Card is corrupt.	Use new data card.
Data Card is not found by software	Data card is blank, no new data.	Data Card is OK, send to patient to download data from SmartLink Module.
Data Card found by SmartLink software but did not download	No new data has been copied from a SmartLink Module.	Data Card is OK, send to patient to download data from SmartLink Module.

## Module

Problem	Probable Cause	Remedy
CPAP does not recognize attached module	The CPAP is not connected to a power source.	Connect the CPAP to an appropriate power source and verify that the module is recognized.
	Module is not fully engaged to the CPAP.	Remove the module from the CPAP, realign the primary data connector with the circular connector of the CPAP, and press firmly until the module is locked into position.
CPAP does not recognize inserted data card when inserted into the module card slot.	Card not oriented properly.	Reference the symbol on the top of the module for proper card orientation.
	Card not fully inserted in the module card slot.	Press the card fully into the card slot until the card clicks into position.
	Module is not fully attached to the CPAP.	Remove the module from the CPAP, realign the primary data connector with the circular connector of the CPAP, and press firmly until the module is locked into position.
Card error on LCD	Wrong card type.	Only use DeVilbiss SmartLink Data Cards. The card must be an SD card with FAT16 format.
	Card removed during data transfer.	Card must not be removed during data transfer. Data transfer interruption may corrupt the card and prevent the card from working correctly. Reinsert card and allow data transfer to continue until completed.
	Card corrupt or faulty.	Return card to provider for replacement.
Settings error on LCD	Card removed during data transfer.	Card must not be removed during data transfer. Data transfer interruption may corrupt the card and prevent the card from working correctly. Reinsert card and allow data transfer to continue until completed. If this error is not corrected, contact provider for resolution.
	Card corrupt or faulty.	Return card to provider for replacement. If this error is not corrected, contact provider for resolution.
Setup required message on LCD	Module not setup properly.	Notify provider of the situation and follow provider instructions.
Oximeter data (SpO2 and PR) not displayed.	Oximeter not properly connected to the SmartLink module.	Verify connection between the oximeter connector and module oximeter connection port.
	Oximeter probe not properly attached to user.	Reference the operating instructions that came with the Oximeter.
CPAP is not fully functional.	Module may be defective.	Remove the module from the CPAP and verify that the CPAP is fully functional. If it functions properly without the module connected, the module may be defective. Contact provider. If it still is not fully functional, consult the CPAP instruction guide for troubleshooting.

## PARTS

DV5M	SmartLink Module
DV5SS	SmartLink Software Kit (CD with software & manual, USB extension cable, SD card, card reader)
DV5FC-10	SmartLink Data Card (10 pack)
DV5FC-1	SmartLink Data Card (single)
DV5CR	SmartLink Card Reader
DV51D-694	USB Extension Cable
DV5M-697	Oximeter Kit (oximeter and 25-use finger probe)
8000AA	Oximeter Finger Clip Sensor
8000Q	Oximeter Ear Clip Sensor
9054D-699	Nonin® nVISION® Software
A-DV5M	SmartLink Module Instruction Manual
A-DV5SS	SmartLink Clinical Software Instruction Manual

## SPECIFICATIONS

### DeVilbiss® SmartLink™ Module

Size ..... 3.9" H x 3.1" W x 1.4" D (9.9cm H x 7.9cm W x 3.6cm D)

Weight ..... 0.30 lbs. (0.14 Kg)

### Electrical Rating

Power Consumption Max (SmartLink module alone)..... 4 Watts

Voltage and Current ..... Power supplied by CPAP

Class II Equipment; Type BF Applied Parts; Continuous Operation

### Operating Conditions

Temperature Range ..... 41°F to 104°F (5°C to 40°C)

Humidity Range ..... 0 to 95% R.H. non-condensing

Atmospheric Pressure Range ..... Sea level to 9,000' (2.743 M)

### Transportation and Storage Conditions

Temperature Range: ..... -40°F to 158°F (-40°C to +70°C)

Humidity Range: ..... 0 to 95% R.H. non-condensing

### DeVilbiss SmartLink™ Desktop Software Specifications

Operating System..... Microsoft® Windows XP™ service pack 2 or higher OR  
Microsoft® Windows Vista™ service pack 1 or higher

Adobe® Reader® ..... 6.0 or greater

RAM..... 512 MB or greater

Hard Drive Free Space..... 150 MB

Hardware: .....SD Memory card reader supporting FAT16 file format OR  
USB port for auxiliary card reader

## DECLARATION OF CONFORMITY

Manufacturer: Sunrise Medical dba DeVilbiss Healthcare  
Address: 100 DeVilbiss Drive  
Somerset, Pennsylvania 15501-2125 USA

Product Designation: CPAP Accessory

Type, Model: DeVilbiss® SmartLink™ Module  
Model DV5M

We herewith declare that the above-mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

Quality System Standards Applied: ISO9001/ISO13485: 2003

Notified Body: TÜV NORD

MDD: Annex II Applied

Safety Standards Applied: IEC 60601-1: 1988 + A1 1995,  
IEC 60601-1-4: 2000,  
ISO 14971: 2000,  
EN ISO 17510: 2000

EMC Compliance to: IEC 60601-1-2: 2001

Authorized Representative: Sunrise Medical Ltd.  
Sunrise Business Park  
High Street  
Wollaston, West Midlands DY8 4PS  
ENGLAND  
44-138-444-6688

®Registered U.S. Patent and Trademark Office and other countries.

Microsoft® Windows® XP™, Windows Vista®, and Internet Explorer™ are registered trademarks of Microsoft Corporation

Adobe® Reader® are registered trademarks of Adobe Systems.

Nonin® nVISION® are registered trademarks of Nonin.

DeVilbiss®, SmartLink™, IntelliPAP™, and SleepCube™ are registered trademarks of DeVilbiss Healthcare.

## LICENSE

One registered copy of the DeVilbiss® SmartLink™ Desktop Software may either be used by a single person who uses the software personally on one computer or installed on a single workstation used non-simultaneously by multiple people.

Full license details are included in the SmartLink Desktop Software install CD in text files named: lgpl-3\_0.txt, LICENSE-2\_0.txt and NOTICE.txt

## LIMITED WARRANTY

### DeVilbiss® SmartLink™ Desktop Software

DeVilbiss warrants that the software installation CD is functional and free of defects. In such case as the user finds the CD to be unusable or defective, it will be replaced free of charge by DeVilbiss

This software is provided "AS IS," without a warranty of any kind. ALL EXPRESS OR IMPLIED CONDITIONS, REPRESENTATIONS AND WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT, ARE HEREBY EXCLUDED. DEVILBISS HEALTHCARE (DEVILBISS) AND ITS LICENSORS SHALL NOT BE LIABLE FOR ANY DAMAGES SUFFERED BY LICENSEE AS A RESULT OF USING, MODIFYING OR DISTRIBUTING THIS SOFTWARE OR ITS DERIVATIVES. IN NO EVENT WILL DEVILBISS OR ITS LICENSORS BE LIABLE FOR ANY LOST REVENUE, PROFIT OR DATA, OR FOR DIRECT, INDIRECT, SPECIAL, CONSEQUENTIAL, INCIDENTAL OR PUNITIVE DAMAGES, HOWEVER CAUSED AND REGARDLESS OF THE THEORY OF LIABILITY, ARISING OUT OF THE USE OF OR INABILITY TO USE THIS SOFTWARE, EVEN IF DEVILBISS HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

### Governing Law

This agreement shall be governed by the laws of the Commonwealth of Pennsylvania, USA.

APPENDIX

Record Setup Screens

My Profile Setup

Required Fields are marked

User Information	
Name	System Administrator
Address	
City	
State/Province	
Postal Code	
Country	
Phone Number	
Email Address	
Contact	
Low Use Threshold	4 Hours
Percentile Pressure	90 Percentile
<div>SaveCancel</div>	

Patient Record Setup

Required Fields are marked

Patient Information	
First Name	
Last Name	
Patient ID	
Address	
City	
State/Province	
Postal Code	
Country	
Phone Number	
Birth Date	
Height	
Weight	
Gender	<input checked="" type="radio"/> Male <input type="radio"/> Female
Email Address	
Timezone	(GMT -05:00) America/New_York
Interface Type/Size	
Interface Accessories	
Payor Information	
Payor Name	
Insurance Id	
Patient Record	
Referring Physician Information	
Name	
Group/Practice	
Address	
City	
State/Province	
Postal Code	
Contact	
<div>SaveCancel</div>	

## Device Record Setup

### Creating a Stand-alone Record

Required Fields are marked

Device Details	
Serial Number	<input type="text"/>
Comments	<input type="text"/>
<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

### Creating a Study Record

Study Details	
Select a Patient	Little Blinken * Archived patients will be restored upon save.
Select a Device	Select a Device
Study Length	60 Days
Comments	<input type="text"/>
<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

## Study Record Setup

### Therapy Settings Tab

#### AutoAdjust models in AutoAdjust mode

Study Details	
Select a Patient	Little Blinken * Archived patients will be restored upon save.
Select a Device	HD001020 (IntelliPAP AutoAdjust)
<b>Therapy Settings</b> <b>Device Settings</b> <b>Respiratory Event Detection</b>	
Therapy Mode	AutoAdjust
Pressure Unit	cmH2O
Upper Pressure Limit	15
Lower Pressure Limit	5
AutoAdjust Delay Time	20 (minutes)
Use Split Night Delay	<input type="checkbox"/>
Split Night Delay Time	240 (minutes)
Number of Nights	1
CPAP Pressure	10
Delay Pressure	4
Delay Time	20 (minutes)
Study Length	60 Days
Comments	<input type="text"/>
<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

#### Standard models and AutoAdjust models in CPAP mode

Study Details	
Patient Name	Sleepy Head
Device	HD001030 (IntelliPAP AutoAdjust)
<b>Device Settings</b>	
<b>Therapy Settings</b> <b>Device Settings</b> <b>Respiratory Event Detection</b>	
Therapy Mode	CPAP
Pressure Unit	cmH2O
Upper Pressure Limit	15
Lower Pressure Limit	5
AutoAdjust Delay Time	20 minutes
Use Split Night Delay	<input type="checkbox"/>
Split Night Delay Time	240 minutes
Number of Nights	1
CPAP Pressure	10
Delay Pressure	4
Delay Time	20 minutes
Study Length	60 Days
Comments	<input type="text"/>
<input type="button" value="Save"/> <input type="button" value="Cancel"/>	

## Study Setup

### Device Settings Tab

#### Standard and AutoAdjust models

### Respiratory Event Detection Tab

#### AutoAdjust models ONLY

## Therapy Setting Options and Defaults

### Therapy Settings tab

#### Standard models and AutoAdjust models in CPAP mode

Setting Option	Default Value
Pressure Unit	cmH2O *
CPAP Pressure	10.0
Delay Pressure	4.0
Delay Time	20
Delay Time Patient Access	Checked

#### AutoAdjust models in AutoAdjust mode

Setting Option	Default Value
Pressure Unit	cmH2O *
Upper Pressure Limit	15.0
Lower Pressure Limit	5.0
AutoAdjust Delay Time	20
Delay Time Patient Access	Checked
Use Split Night Delay	Unchecked
Split Night Delay (only available if Use Split Night Delay is Checked)	240
Number of Nights (only available if Use Split Night Delay is Checked)	1
Study Length	7

### Device Settings tab:

#### Standard and AutoAdjust models

Setting Options	Default Value
LCD Language	English *
Tubing Length	1.8m (6ft)
Patient Access	Checked
Enable Auto ON	Checked
Enable Auto OFF	Unchecked
Enable Notifications	Checked
Enable Mask Fit Check	Checked

### Respiratory Event Detection tab:

#### AutoAdjust models ONLY


Setting Option	Default Value
Apnea Flow Reduction	10
Apnea Duration	10
Hypopnea Flow Reduction	50
Hypopnea Duration	10

\* DV5X-NE Series SleepCube devices have factory defaults of hPa for Pressure Unit and German for LCD Language.

## Sample Reports


## DV54 AutoAdjust

## Summary Report Summary view- total hours of usage

DV54 AutoAdjust

Summary Report Detail view-ON/OFF usage

## DV54 AutoAdjust



## One-Night Report

Study Details			
Patient Name	Little Winken	Device Type	HD002853 (IntelliPAP AutoAdjust)
Start Date	6/28/08	Study Length	7

Summary 1-Night

Sessions: &lt;&lt; 1: 6/28/08 10:47 PM - 6/29/08 6:41 AM &gt;&gt; Refresh Export to PDF

6/28/08 10:47 PM - 6/29/08 6:41 AM

## Events Summary (Delay Mode)

Time: 0.0 Hours  
AHI: -  
AI: -  
HI: -  
NI: -  
Snores: -  
Mask Leak Time: -  
Exhale Puffs: -  
Average Pressure: -  
Average Leak: -  
Avg. Est. Tidal Vol.: -  
Avg. Breath Rate: -


## Events Summary (AutoAdjust Mode)

Time: 8 Hours  
AHI: 5 Per Hour  
AI: 0 Per Hour  
HI: 5 Per Hour  
NI: 0.4 Per Hour  
Snores: 93 Events  
Mask Leak Time: 0 % of total  
Exhale Puffs: 198.2 Per Hour  
Average Pressure: 7.4 cmH2O  
Average Leak: 32.3 L/min  
Avg. Est. Tidal Vol.: 534.9 mL  
Avg. Breath Rate: 18.4 bpm

## DV5I Standard

### Summary Report

## Patient Notes

### System Notes (Auto-Entry)

All Patient Notes include the date of entry and source of entry.

Patient Note—System Note	Event Trigger
First data received for study	The first successful transfer of therapy data from the CPAP to the software occurred.
SmartLink Data Card has been programmed with updated settings. Verification of settings is pending.	New therapy settings were transferred from a study record to a data card and the software is waiting for confirmation from the CPAP that the changes were completed.
Settings update occurred on <date from SD file>, and validated on <today's date>	The CPAP successfully completed the therapy setting changes sent by the software via a data card.
Device <serial number> was assigned to patient	A device record was assigned to a patient record, renewing or creating a study record.
Device <serial number> was un-assigned from patient	A device record was removed from a patient's study record.
Errors in data import. Some data may have been lost.	A transfer attempt from the module or software to the data card did not complete successfully.
New data has been uploaded	A transfer attempt from the data card to the module or software completed successfully.
Download contains one or more faulty records	A transfer attempt was not entirely successful.

## SmartLink Messages - shown on CPAP LCD Display

Oximeter Connected	Oximeter is connected and logging has started. This displays for five seconds; then switches to the pulse oximeter menu (SpO2/PR).
Oximeter Removed	Oximeter is removed. This displays intermittently for one minute.
SmartLink Module Ready for Use	Module is connected correctly. Displays for ten seconds.
SmartLink Menu Press ◀▶ To Exit	Access and exit the menu by pressing the LEFT and RIGHT arrow keys at the same time.
Quick View Menu Press ▼	Access the Quick View menu by pressing the DOWN arrow key.
Average Pressure 30 day avg: 10.1	The average pressure provided during the selected time period is displayed for 1, 7, or 30 days.
AHI 30 day avg: 2.0	The average AHI (Apnea/Hypopnea Index: the number of apneas and hypopneas detected per hour) is displayed for 1, 7, or 30 days.
Leak % 30 day avg: 5.0	The average leak percentage (percentage of therapy time that the system leak was above 95 LPM) is displayed for 1, 7, or 30 days.
NRI 30 day avg: 1.0	The average NRI (Non-Responding Event Index: number of non-responding events-exhale puff, leak, and mixed events-detected per hour) is displayed for 1, 7, or 30 days.
EPI 30 day avg: 1.0	The average EPI (Exhale Puff Index: number of exhale puff events detected per hour) is displayed for 1, 7, or 30 days.
Hours Per Day 30 day avg: 8.0	The average number of hours of therapy provided during the past 1, 7, or 30 days is displayed.
Adherence % 100% (1234/1234)	The percentage of days in which the system was used for at least 4 hours per day is displayed for 1, 7, or 30 days.
Exit Menu Press ▲	Press the UP arrow key to exit the Quick View menu.
Volume 100% (Adjust ▲▼)	Volume level for audible alert. Press the UP or DOWN arrow key to adjust.
Serial Number CM123456	Serial number of the module.
Firmware Version V0.06 25/04/2008	Version of operating firmware in the SmartLink module.
SpO2: 100 PR: 101	Pulse oximetry (SpO2) and pulse rate (PR) are displayed if a pulse oximeter is used with the system. Values read *** if a pulse oximeter is not connected.
Card Detected Please Wait . . .	Data card is inserted and data is being transferred. This displays until the transfer is complete.
Card Detected Ready For Use	Data card is inserted and there is no new data to transfer. This displays for ten seconds.
Card Transfer Successful	Data transfer to data card complete. This displays for ten seconds.
Card Error Contact Provider	Error occurred during data transfer. This displays until acknowledged or for 15 minutes.

Setting Error  
Contact Provider

Error occurred while updating device settings. This displays until acknowledged or for 15 minutes.

Card Removed  
Send to Provider

Data card is removed from the module. This displays for ten seconds.

Setup Required  
Contact Provider

Module and CPAP are not matched. The module is unable to record data; however, the CPAP will continue to provide therapy. Displays until the module is removed or data is cleared by holding Down Arrow key for 10 seconds.

E00 Module Fault  
Contact Provider

A module fault has occurred (Clock Error, Memory Error, Setting Error, etc).

Firmware Update  
Please Wait . . .

Firmware is being updated. This displays until update is complete and system operation is normal.

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The screenshot shows the 'Detalles de dispositivo' form. It has a header 'Los campos obligatorios aparecen indicados'. Below it, there are two main sections: 'Número de serie' with a text input field, and 'Comentarios' with a text area. At the bottom, there are two buttons: 'Guardar' and 'Cancelar'.

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The screenshot shows the 'Detalles de dispositivo' form. The 'Número de serie' field is filled with 'ID002580'. The 'Comentarios' field is empty. At the bottom, there are two buttons: 'Guardar' and 'Cancelar'.

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The screenshot shows a navigation menu with 'Appareils' selected. Under 'Appareils', there are three sub-items: 'Tous les appareils', 'Appareils disponibles', and 'Appareils affectés'. The 'Nouvel appareil' button is located at the bottom of the 'Appareils' section. A dropdown menu for '10 éléments par page' is also visible.

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The screenshot shows the 'Détails de l'appareil' form with the 'Numéro de série' field filled with 'ID002580'. The 'Observations' text area is empty. At the bottom, there are 'Enregistrer' and 'Annuler' buttons.

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The screenshot shows the 'Dettagli dispositivo' form. It has a header 'I campi richiesti sono contrassegnati'. Below the header, there are two main sections: 'Numero di serie' with a text input field, and 'Comments (Commenti)' with a text area. At the bottom, there are two buttons: 'Salva' and 'Annulla'.

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The screenshot shows the 'Dettagli dispositivo' form. It has a header 'I campi richiesti sono contrassegnati'. Below the header, there are two main sections: 'Numero di serie' with a text input field containing 'ID002580', and 'Comments (Commenti)' with a text area. At the bottom, there are two buttons: 'Salva' and 'Annulla'.

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®Registrado en la Oficina de patentes y marcas (Patent and Trademark Office) estadounidense y en otros países.
®Eingetragen beim U.S. Patent and Trademark Office und in anderen Ländern.
**Sunrise Medical Pty. Limited**